

Autoclaves and Accessories			
Product Code	Product Name	Description	List Price
T50.000.00	SX-500	Autoclave Single phase, AC120V, 50/60Hz, 58 L (Effective Volume 50 L), Max. Temp 135° C, includes two ACA-315B short baskets (accommodates three)	\$13,373.35
T70.000.00	SX-700	Autoclave Single phase, AC230V(207V~240V), 50/60Hz, 69 liters (60L effective capacity), warming and heating cycle, two baskets ACA-700B included	\$16,225.88
T01.000.00	ACA-315T	STAINLESS TALL BASKET for SX-500 / 300 dia. X557 ht.mm /	\$596.00
T01.002.00	ACA-315K	STAINLESS BUCKET / 299 dia.X279 ht.mm for SX-500 (up to 2pcs installable)	\$916.34
T01.003.00	ACA-315B	STAINLESS BASKET/ 300 dia.X182 ht.mm for SX-500 (up 3pcs. installable)	\$571.39
T01.018.00	ACA-325TP	STAINLESS LONG BASKET for SX-500 (300dia. X491 ht.mm)	\$915.00
T01.010.00	ACA-700T	STAINLESS TALL BASKET for SX-700 / 345dia.X561 ht.mm	\$1,056.00
T01.012.00	ACA-700K	STAINLESS BUCKET / 344 dia.X279 ht.mm for SX-700	\$999.42
T01.009.00	ACA-700B	ACA-700B-STAINLESS BASKET / 345 dia.X181 ht. mm for SX-700	\$687.69
T01.020.00	ACA-710TP	STAINLESS LONG BASKET for SX-700 345dia.X527 ht.mm	\$1,374.00
T02.015.00	Lid Gasket for SX-500	Lid Gasket of SX-300/SX-500	\$849.00
T02.016.00	Lid Gasket for SX-700	Lid Gasket for SX-700	\$963.00
T02.004.00	Printer for SX-500/700	PRINTER for SX-500/700 T	\$2,077.59
T02.011.01	External temperature sensor protector kit SX-500	Protector kit for external temperature sensor for articles to be sterilized for SX-500 autoclave	\$642.09
T02.011.02	External temperature sensor protector kit SX-700	Protector kit for external temperature sensor for articles to be sterilized for SX-700 autoclave	\$642.09
T02.109.00	EXTERNAL SENSOR for SX-500/700	External temperature sensor for SX-Series Autoclave for sterilization media (part only, labor sold separate)	\$839.24

Centrifuges and Accessories			
Product Code	Product Name	Description	List Price
T31.000.00	MDX-310	High Speed Refrigerated Centrifuge; 15000 rpm; 21,890G	\$12,951.30
T31.100.00	CA-300	Rotor rack bucket required for MDX-310	\$3,132.60
T01.501.00	AR015-24	Rotor Rack 1.5ml/2.0mlx24/15,000rpm	\$1,023.50
T01.509.00	AR510-04 (conical)	Rotor Rack 15mlTC/50mlTCx4/15,000rpm	\$1,927.40
T01.503.00	AR005-24	ROTOR RACK / 0.5ml x 24 / 15,000rpm	\$1,056.32
T01.502.00	AR002-64	ROTOR RACK / PCR 0.2ml x 8tubes x 8 strips / 15,000rpm	\$1,084.11
T01.510.00	AR501-04	ROTOR RACK / 14ml/50mlx4/15,000rpm	\$2,167.75
T30.928.00	AR150-08 (conical)	ROTOR RACK / 15mlTC x 8 / 15,000rpm for TMA-300	\$2,259.16
T01.511.00	PCR96-02	ROTOR RACK / Adapter A96-01PC required	\$3,313.04
T01.512.00	A96-01PC	96-well plate adapter working only for non-skirt types / PCR Plate x 2/14,000rpm (Qty.2)	\$2,108.76
T01.507.00	AR015-SC24 (spin column)	ROTOR RACK for Spin Columns / 1.5ml/2.0mlx24/15,000rpm	\$1,505.35
T01.504.00	AR050-12 (conical)	Rotor Rack 5mlx12/15,000rpm	\$2,101.83
T01.505.00	AR141-08	Rotor Rack 14ml/8/15,000rpm	\$1,717.91
T01.506.00	AR300-06	ROTOR RACK / 30ml x 6 / 15,000rpm	\$1,657.19
T30.910.00	CA-1	Fixed Angle Rotor, 1.5 / 2.0 mL x 18	\$3,502.62
T30.915.00	CA-10	Fixed Angle Rotor; 2 mL x 24	\$3,130.89
T30.941.00	CA-15	FIXED ANGLE ROTOR 2.0mlX18, 0.5&0.2X12/15,000rpm	\$4,712.04
T30.912.00	CA-5	Fixed Angle Rotor; 2 mL x 24	\$4,335.08
T30.929.00	CA-6	FIXED ANGLE ROTOR mlX36/15,000rpm	\$6,473.82
T30.940.00	CA-12	FIXED ANGLE ROTOR mlX12/15,000rpm	\$6,141.36
T30.927.00	CA-17	FIXED ANGLE ROTOR 50ml or 15ml x 4/15,000rpm for MDX-310	\$9,758.21
T30.200.00	MultiSpin	TOMY MultiSpin Cordless Battery-Powered Mini Centrifuge; 6,000RPM x 8 1.5 / 2.0 mL Microtubes (2,960G)	\$390.00

Autoclaves Extra Services			
Product Code	Product Name	Description	List Price
T99.002.00	SX-Series Altitude Factory Reprogramming	Required above 1,000 meters. Please specify altitude.	\$1,057.56
T99.300.00	3 Year Extended Warranty for Tomy Autoclave (SX 700 series)	Extension of standard 1 year warranty to 3 years total (extension of 2 additional years)	\$699.00
T99.300.01	3 Year Extended Warranty for Tomy Autoclave (SX 500 series)	Extension of standard 1 year warranty to 3 years total (extension of 2 additional years)	\$597.00
T99.500.00	5 Year Extended Warranty for Tomy Autoclave (SX-700 series)	Extension of standard 1 year warranty to 5 years total (extension of 4 additional years)	\$1,398.00
T99.500.01	5 Year Extended Warranty for Tomy Autoclave (SX-500 series)	Extension of standard 1 year warranty to 5 years total (extension of 4 additional years)	\$1,194.00
T11.201.92	TOMY - Palletizing	Wrapping, Strapping, and Palletizing TOMY Autoclave for safe shipment. Required for SX-500 and SX-700	\$106.00