

Improve sterilization workflow, usability, and efficiency with the following accessories.



Tall Basket

ACA-315T (SX-500) / ACA-700T (SX-700)

- Combination perforated steel with a solid-bottom basket that takes up the entirety of the chamber.
- For autoclaving large items such as bioreactors/fermenters and bio waste.
- Doubles as wastebasket when not in use (use autoclave bag).
- Autoclave accommodates 1.
- Dimensions (Diameter x Height):
300 x 557 mm (ACA-315T)
345 x 561 mm (ACA-700T)



Stainless Buckets

ACA-315K (SX-500) / ACA-700K (SX-700)

- Stackable solid stainless steel bucket for autoclaving liquids and medium-size items.
- Prevents liquids from spilling into autoclave chamber and damaging components
- Autoclave accommodates 2.
- Dimensions (Diameter x Height):
299 x 279 mm (ACA-315K)
344 x 279 mm (ACA-700K)



Short Basket

ACA-315B (SX-500) / ACA-700B (SX-700)

- Autoclave includes 2 baskets; add a third basket to increase loading capacity by 50%.
- For small glassware and instruments (liquids not recommended).
- Stackable perforated steel basket to allow for easy steam penetration.
- Dimensions (Diameter x Height):
300 x 182 mm (ACA-315B)
345 x 181 mm (ACA-700B)



Replacement Lid Gasket

For SX-500 and SX-700

- For replacement every 1,000 operating hours.
- Helps to maintain proper function of autoclave.
- Prevents steam leakage into autoclave and damage to electronics, other sensitive equipment (that could lead to thousands of dollars in repair damage).

We can assist with questions about fit, capacity and compatibility-
Please [contact us!](#)

(858) 800-3900 | tomy@amuzainc.com | tomy.amuzainc.com

